

Independent Living. Dementia Care. Lone Working.

Indicare: The Wearable Personal Safety Wristwatch

A stylish wristwatch for Elder Care, Independent Living, Dementia & Lone Working.





- Automatic fall detection
- SOS Button
- Mobile phone in an analog wristwatch
- Beacon sensitive
- Locates Users both outdoors & indoors!
- Long battery life
- Automatic location and battery alarms
- Out of bounds alerts
- No movement alarm
- Tailormade configurations
- Embedded GPS

The Halo Indicare wristwatch forms part of our tailorable service for the safety of people. Indicare is perfect for everyone from the Independent Living seeking peace-of mind, say, in case of a fall, to dementia or Alzheimer's sufferers. Its rugged design makes it ideal for Lone Workers too.

It's your call! Indicare does as much or as little as you need it to:

If you fall: Indicare lets our monitoring station know so we can get you help, fast

If you need phone a friend: Indicare lets you talk to up to 7 contacts on its inbuilt loud speakerphone

If you need urgent help: Press & hold the crown button to send an SOS to our Monitoring Centre

Contact the wearer: Up to 7 contacts may phone Indicare and chat with the wearer

Indoor positioning: Optional beacons allow Indicare to know where it is indoors

GPS positioning: GPS function allows the User to be located globally in an emergency

Eldercare -Activity monitoring: Flags changes in activity which may be an indication of failing health

Wander prevention: Be alerted immediately if a dementia sufferer leaves a given location

Indicare is a highly versatile system with too many variables to list. Therefore, we highly recommend discussing your needs with a Halo Group sales professional who can tailor a solution to fit your precise needs.

308 Harbour Point Business Park Little Island Cork T45 EW73 Ireland.



Independent Living. Dementia Care. Lone Working.

Indicare: The Wearable Personal Safety Wristwatch

Key features are:

- Aestethically pleasing wristwatch with normal hour and minute hands
- 3D-accelerometer, support for wireless sensors, on-board GPS, QUAD band cellular processor

Industry benchmark for footprint and functionality:

Mobile comms	Modem type	Quad band GSM/GPRS/EDGE
	Communication	SMS, FTP, HTTP, TCP, UDP, USSD
GNSS	Receiver type	GPS/GLONASS receiver, SBAS
		WAAS/EGNOS/MSAS support
	Sensitivity	Acquisition -148 dBm, navigation -162 dBm
	A-GPS	A-GPS assistance (up to 28 days)
Short range comms	BLE	Short range comms BLE Bluetooth Low Energy
		communication with home and mini beacons
мси	Processor	ARM Cortex-M3, 512 kbytes (total code
		memory)
	RAM & flash	64 kbytes (total RAM), 8 Mbytes (total data
		flash)
Audio	Microphone / speaker	Adjustable gain and volume control
User interface	Wearer profiles	Tracking (dementia), Auto, Easy and Full
		profiles
	Adjustable UI	Selected profile enables/disables UI features
	Watch	Analog. Two discretely software controllable
		hands
	Buttons	Crown, emergency call, speed dial and
		yes/no/up buttons
	LEDs	3 user interface / status LEDs
	Other	Two row OLED display, speaker, vibrator
Charging		Easy to use integrated charging dock in the
		home beacon
Sensors	Onboard	3D accelerometer, temperature
Battery		Li-ion 500 mAh
Operation time	GSM always on	Typical 2 - 7 days, depending on use profile and
		conditions
	GSM on when requ <mark>ired</mark>	Up to 50 days, depending on use profile and
		conditions
Software features	Data communicatio <mark>n</mark>	SMS, USSD and GPRS with message buffering
		on data flash
	Over-the-air update <mark>s</mark>	Configuration, parameters, safe zone, firmware
		and boot loader
	Safe zone functiona <mark>lity</mark>	Polygon or circle safe zones with validity time
		attributes
	User initiated phon <mark>e calls</mark>	To alarm call centre and up to 7 speed dial
		numbers
	Receiving phone cal <mark>ls</mark>	Automatically or by user accepted calls, white
		list access control
	Location alerts	Automatic location alerts for dementia person
		leaving home beacon

308 Harbour Point Business Park Little Island Cork T45 EW73 Ireland.





Independent Living. Dementia Care. Lone Working.

		macpendent in Bright and a care in in
	Alarms	Home beacon out of range, safe zone breach,
		battery, no movement, Man down
	Activity statistics	Wearer statistics including duration of very
		active, active, little movement, at home by
		area, out of home periods
	Event reports	Periodic status, help request call, home
		beacon, safe zone, activity statistics, battery
		status and charging reports, location reports
Temperature range	Operation	Use: -20°C - +60°C, S1 charging 0°C - +45°C
Certifications		RoHS, CE, 99/5/EC
Enclosure	Indicare SX1	Water protected (brief submersion), standard
		strap, 22 mm interface
Dimensions	Indicare SX1	Diameter 48 mm, thickness 17 mm, weight 65
		gr. (normal strap)
Home beacon	Charging	Easy to use charging dock
	Mini beacons	Provide location and statistics information
		indoors